
- EN CERVICAL COLLAR
- DE ZERVICALDRITHESE
- FR COLLIER CERVICAL
- ES CUELLO CERVICAL
- IT COLLARE CERVICALE
- NL NEKBRACE
- DK HALSHVÆRSELKRAVE
- FI KALLATIUKKI
- SE MACKKRAGE
- NO CERVICAL KRAGE
- PL SZYJANY KOLANIEZ ORTOPEDYCZNY
- GR ΑΥΞΗΝΙΚΟ ΚΟΛΑΠΟ

ENGLISH

Product Description
A lightweight two-piece collar in polyester, with rigid occipital and mandibular posts.

Usage Information (See diagram)
1. After selecting correct size, snugly back half of collar to back of patient's neck. Height – measure bottom of chin to top of sternum.
2. Apply front half of collar with chin secured in recess. Centre collar to ensure medial alignment. Front half overlaps back half.
3. With hoodlock fasteners, tighten collar with arrow pointing up.

Indications for Use
• Cervical injuries for trauma patients.
• Immobilization for pre-collapse surgery.
• Immobilization conditions which require gross immobilization of the mid cervical spine.

Contraindications
• Severe neck injuries such as cervical vertebral fractures.
• Cervical dislocation with head angulation.
• Soft tissue swelling in the neck.

Warnings and Precautions
We recommend the initial fitting of this brace be conducted by a suitably qualified healthcare professional who will advise the period of use.
• Carefully read all instructions and warnings prior to use.
• Follow all instructions to ensure proper performance of the brace.
• Do not use if fragments, ornaments, girths, creams or any other foreign matter is present on the skin.
• Do not reuse for another patient, due to cross infection and cross contamination.
• Perform regular skin and circulation checks, especially patients with diabetes, vascular deficiencies and neurological conditions.
• Should any adverse reactions occur, please cease use and contact your healthcare professional or provider of this product.
• The elasticity of the brace may be compromised by certain and/or fasteners. To prevent this the hook and loop fasteners should always be fastened where the brace is not being worn or when being washed.

Washing and Care Instructions
1. Close over all hoodlock straps.
2. Clean the front and back panels of the collar with a cloth and mild detergent.
3. Do not iron.
4. Do not stain or tumble dry.
5. Do not use fabric conditioner at any time.
6. Check for any degradation that may affect use.
7. Periodically check hoodlock fasteners are operating securely.

Product Family Composition
Nylon, Cotton, Polypropylene, Polyester

Storage and Transport Conditions
Store in a cool, dry place out of direct sunlight. In the original packaging.

Recycling and Disposal
Packaging and condition parts should be recycled or safely disposed in accordance with local or national laws.

Setious Incident
Report any serious incident to the manufacturer and the competent authority of the EU Member State, or country in which you reside.

Recycling and Disposal
Packaging and condition parts should be recycled or safely disposed in accordance with local or national laws.

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©2021 Promedics Orthopaedics Ltd. All Rights Reserved.
Promedics Orthopaedics Ltd.